

BASSETT HEALTHCARE NETWORK

Date and Time of Verification _____

CSF Verification Form

MR# _____

CLP-Contact provider to confirm testing and collection sequence- Tell provider "You may also receive a call from a Technologist if there is a question about integrity of the sample."

Patient Name _____

Deliver specimens to appropriate departments announcing a STAT CSF

Ordering Provider _____

Verifying Provider _____

Verified By _____

Is there a *14-3-3 Protein (CJD) test Ordered? If yes **STOP** all Testing Must be sent to ARUP lab.

Have all test been released? YES NO

Verify Collection Sequence: (What order were tubes collected) _____

Test	Tube Number(s) testing to be performed on
*14-3-3 Protein (CJD)	If ordered STOP all Testing Must be sent to ARUP lab
Cell Count (Heme)	
Total Protein (Chem)	
Glucose (Chem)	
Lactate (Chem)	
Cultures/other Micro Testing (Micro)	
Cryptococcal Antigen (Micro)	
Cytology (Cyto)	
Herpes Simplex Virus (HSV) (Send-out)	
Varicella-Coster Virus (VZV) (Send-out)	
VDRL (Send-out)	
Viral Culture (Send-out)	
Lyme (Send-out)	
Freeze and Hold (CLP)	
Other-Specify Testing	

<u>Neurology CSF Panel</u>	<u>Default Tube number</u>
CSF Cell count with differential (LAB212)	1 & 4
Fluid Culture/Gram Stain (LAB269)	2
AFB culture and stain (LAB877)	2
Glucose CSF (LAB185)	3
Protein CSF (LAB195)	3

<u>Standard CSF Panel</u>	<u>Default Tube number</u>
CSF Cell count with differential (LAB212)	3
Glucose CSF (LAB185)	1
Protein CSF (LAB195)	1
Fluid Culture/Gram Stain (LAB269)	2
Cytology, Non Gyn (LAB13)	4

Note: Providers have the option to change the tests ordered on each tube.

- Cerebrospinal fluid must be submitted in the proper collection tubes provided with the Lumbar Puncture Tray (warehouse item #0000115).
- ASCP recommends the following testing by order of tubes:
 - **Tube #1** will be used for chemistries and immunology
 - **Tube #2** will be used for microbiology and cultures
 - **Tube #3** will be used for cell counts and differentials
 - **Tube #4** if collected, will be used for cytology

CONSULT LABORATORY MANUAL ON BASSETT INTRANET FOR PROPER COLLECTION PROCEDURES.